



PEDIATRIC ZIRCONIA CROWNS





New Generation **Pediatric Zirconia** Crowns

“One Session So Much Happiness”

www.profzrcrown.com • 0850 241 90 36 •  [prof.zrcrown](https://www.instagram.com/prof.zrcrown)



Why ProfZrCrown?

As a result of ProfTeknoloji's hundreds of model scans, 2 years of R&D work and more than 300 case studies ProfZrCrowns have been produced.

- ✓ Innovative Spacious and Inclusive Interi-
- ✓ Aesthetic Design
- ✓ Increased Stamina
- ✓ Ease of Application
- ✓ Size Variety
- ✓ Innovative Turkish Design









ProfZrCrown offers a new coding technique that does not disturb the language. Specific with the numbering system and inscription in the crown, easier tooth selection and children's language habit in the embossed part of the post-application text eliminates the possibility.

Case: Dr. Ozan Baris AYDIN



Sizing Information

ProfZrCrowns offer 5 different size options in your applications and are available in 3 sizes in our Small, Medium, Large standard boxes. It offers an easier application option with our XSmall and XLarge options when needed in your different cases.

Central			Lateral			Max. Canine			Man. Canine		
											
1-	1-	XSmall	2-	2-		3-	3-	XSmall	3-	3-	
1S	1S	Small	2S	2S		3S	3S	Small	3S	3S	
1M	1M	Medium	2M	2M		3M	3M	Medium	3M	3M	
1L	1L	Large	2L	2L		3L	3L	Large	3L	3L	
1+	1+	XLarge	2+	2+		3+	3+	XLarge	3+	3+	

SIZING GUIDE

Central



	X _{mm}	Y _{mm}	
xS	5.60	6.00	
S	5.75	6.20	
M	6.02	6.50	
L	6.35	6.85	
xL	6.70	7.25	

Lateral




	X _{mm}	Y _{mm}	
xS	5.05	5.20	
S	5.25	5.32	
M	5.63	5.60	
L	5.94	5.90	
xL	6.26	6.25	

Max. Canine



	X _{mm}	Y _{mm}	
xS	5.82	6.40	
S	6.05	6.55	
M	6.53	6.88	
L	6.60	7.16	
xL	6.94	7.50	

Man. Canine



	X _{mm}	Y _{mm}	
xS	5.36	5.30	
S	5.58	5.50	
M	5.85	5.70	
L	6.11	6.05	
xL	6.37	6.31	

” PRODUCT DIVERSITY



Online Order



Have different sizing **perfect aesthetics**
and harmony with crowns..



Prof Teknoloji



PRØF ZrCrown



0850 241 90 36

ProfZrCrown children's zirconium crowns completely in ProfTechnology production facilities is produced.



CONTACT

 0850 241 90 36

 0552 819 69 76

 www.profzrcrown.com

 info@profzrcrown.com